RANOX
FOOD DIAGNOSTICS

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
Recovery test of Pork samples

The Recovery values shown in Graph #2 show that when compared to chemical standards the matrix fortified samples have significant recovery values. The matrix and chemical samples were both calculated against calibration curves made with the Randox standards. The recovery was calculated by the values obtained from the matrix spikes divided by the values obtained from the chemical standards.

Graph #2. Avg % Recovery values

- 1 Sulfadiazine
- 2 Sulfadimethoxine
- 3 Sulfaquinoxaline
- 4 Sulfamethazine
- 5 Sulfathiazole
- 6 Sulfisoxazole
- 7 Sulfapyridine
- 8 Sulfamerazine
- 9 Sulfamethoxypyridazine
- 10 Sulfachloropyridazine
- 11 Sulfamethazole
- 12 Sulfadoxine

Evaluation of Sulfonamide Test kit from Randox in Laboratory Study; Centre for Veterinary Drug Residue, Canadian Food Inspection Agency

Jana Kinar*, Johanna Chrusch, Joe Boison, 116 Veterinary Road, Saskatoon SK S7N 2R3

The detection capabilities CCB β for 12 validated antibiotics (enrofloxacin, difloxacin, ceftiofur, Desfuroyl Ceftiofur Cysteine Disulfide, thiamphenicol, florfenicol, tylosin, tilmicosin, streptomycin, dihydrostreptomycin, tetracycline, doxycycline) were all lower than the respective MRLs in muscle from different animal origins (bovine, ovine, porcine, poultry).

Evaluation and validation of a multi-residue method based on biochip technology for the simultaneous screening of six families of antibiotics in muscle and aquaculture products; Valérie Gaudin, Celine Hedou, Christophe Soumet & Eric Verdon

<http://dx.doi.org/10.1080/19440049.2015.1125529>

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